

# Otyker orthopaedics

Otyker Orthopaedics Technology Co. Ltd.

Phone: +852 8191 3424

Email: [contact@otyker.co](mailto:contact@otyker.co)

Skype: otykerco

Address: 212Rm, No.592 Tong He Road,  
Qingdao 266035 China



## LARGE FRAGMENT LOCKING SCREW RACK

Large Fragment Locking Screw Rack. Very nice quality.

Screw Length: 10-60mm/10-50mm/30-60mm  
Locking Screw Length: 10-60mm/16-60mm



## LARGE FRAGMENT SCREW RACK 4.5/6.5MM

Excellent quality - Large Fragment Screw Rack

4.5mm Screw Length: 20-60mm  
6.5mm Screw Length: 45-80mm



## LARGE FRAGMENT SCREW RACK-II

High quality Large Fragment Screw Rack-II

4.5mm Screw Length: 18-66mm  
6.5mm Screw Length: 30-95mm



## LARGE FRAGMENT SCREW RACK-III

Large Fragment Screw Rack-III. Material made of PPSU.

4.5mm Screw Length: 20-60mm  
6.5mm Screw Length: 35-90mm



## ORTHOPAEDIC INSTRUMENT CASE WITH 3.5/4.0MM SCREW RACK

Orthopaedic screw/instrument case - 3.5/40mm screw rack aluminium. Made of excellent quality aluminium - very strong. The case has two trays. Each tray has rubbers instrument holder that can be cut by your needs to hold the instruments. Has excellent customers review. Features of case with 3.5/4.0 screw rack are listed below.

Dimensions: 15" x 9" x 4" (L x W x H)  
Screw rack removable, with holes market:  
3.5mm 12,14,16,.....42 each size 6 holes  
4.0mm 12,14,16,.....50 each size 6 holes





## ORTHOPAEDIC INSTRUMENT CASE WITH 3.5/4.0MM SCREWS RACK

Excellent quality! Orthopaedic instrument case and boxes with 3.5/4.0mm screw rack. Nice size. The metal case is aluminium and it comes with a very strong cover. Below are features of orthopaedic screw case with rack 3.5/4.0mm.

Dimension of the case is 15.5" x 9.5" x 4" (L x W x H)

Holes for 3.5mm (12,14,16.....42mm) 6 holes each size (96 holes)

Holes for 4.0 mm (12,14,16.....42mm) 6 holes each size (96 holes)

Total holes (3.5 and 4.0) = 192



## ORTHOPAEDIC SCREW CASE WITH RACK 3.5/4.0MM

Veterinary animal orthopaedic instrument case and boxes with screw rack 3.5/4.0mm. It is small in size. Both parts of the case are made from one piece of aluminium, not welded and without screws which makes this case excellent for work.

Dimension of the case (L x W x H) is 10" x 4.5" x 2-3/8"

Dimensions of the screw rack (L x W x H) is 9-3/8" x 4-1/8" x 2-1/3"

Holes for 3.5mm (12,14....40mm) 5 holes each size (120 holes)

Holes for 4.0 mm (96 holes)

Total holes (3.5 and 4.0) = 216



## SCREW RACK 1.5/2.0/2.7MM

Good quality! The screw rack is small in size. Below are the amount of holes it provides for screw size of 1.5mm, 2.0mm and 2.7mm.

Holes for 1.5mm (6, 8,10.....20mm) 5 holes each size (40 holes)

Holes for 2.0 mm (8, 10, 12.....28mm) 6 holes each size (66 holes)

Holes for 2.7mm (10, 12, 14.....30mm) 6 holes each size (66 holes)

Total holes (1.5, 2.0 and 2.7) = 172



## SCREW RACK 3.5/4.0MM

Small Fragment Screw Rack

3.5mm Screw Length: 12-40mm

4.0mm Screw Length: 14-45mm



## SCREW RACK FOR SPINAL PEDICLE SCREW, ROD AND CROSSLINK

Screw Rack for Spinal Pedicle Screw, Rod and Crosslink

Mono Screw, Mono Reduction Screw, Poly Screw, Poly Reduction Screw, Crosslink, Rod



## SMALL FRAGMENT LOCKING SCREW RACK

Small Fragment Locking Screw Rack

Screw Length: 10-60mm/10-50mm/30-60mm

Locking Screw Length: 10-60mm/16-60mm



## SMALL FRAGMENT SCREW RACK 3.5/4.0MM

Small Fragment Screw Rack-II. Very nice size and color.

3.5mm Screw Length: 12-40mm

4.0mm Screw Length: 14-45mm



## SMALL STERILIZATION CASE

Excellent quality! Very good for small autoclaves. It is small in size. It comes in a nice case.

Made of aluminium. Wall is 2mm thick -very strong. It provides two extra trays and rubber instrument holders. Nice blue color lid.

The dimensions of the case is 12" x 6 1/2" x 4" (L x W x H).