

Otyker orthopaedics

Otyker Orthopaedics Technology Co. Ltd.

Phone: +852 8191 3424

Email: contact@otyker.co

Skype: otykerco

Address: 212Rm, No.592 Tong He Road,
Qingdao 266035 China

COMPLETE KIT 1.5/2.0/2.7/3.5/4.0MM



from
\$2999

This is our most popular system! It contains of all necessary standard tools need for orthopedic surgery. It is pack with full range of the instruments and implants. All instruments are high quality. All screws are self tapping-hex head. Made of stainless steel. This kit can be only for veterinary use in the USA.

This kit can be adjusted according to customer's requirements. Instruments and implants included in this complete kit are listed below.

Drill bit 1.1 2 Qty
Drill bit 1.5 2 Qty
Drill bit 2.0 2 Qty
Drill bit 2.5 2 Qty
Drill bit 3.2 2 Qty
Neutral & Loaded Drill Guide - 1.1mm for 1.5mm Screws
Neutral & Loaded Drill Guide - 1.5mm for 2.0mm Screws
Neutral & Loaded Drill Guide - 2.0mm for 2.7mm Screws
Neutral & Loaded Drill Guide - 2.5mm for 3.5mm Screws
Neutral & Loaded Drill Guide - 3.2mm for 4.5mm Screws
2.5mm & 3.5mm Tap Sleeve
3.2mm & 4.0mm Tap Sleeve
2.0mm & 2.7mm Tap sleeve
HA1.5 Tap
HA2.0 Tap
HA2.7 Tap
HA3.5 Tap
HB 4.0 Tap
Depth Gauge
Hex Head Screw Driver SW2.5
Hex head Screwdriver 1.5mm
Plate benders pair-modern design large
Plate benders small
Periosteal Elevator(Round Head)
Broken screw exterioriser
Counter sink large
Counter sink small
Periosteal Elevator (Flat Head)
Reduction Forceps with point
Reduction Forceps
Self-centering forceps 2 Qty
Case aluminum
Screw rack-3.5/4.0mm
Screw rack-1.5/2.0/2.7
Bone plate 1.5mm 5 Qty
Bone plate 2.0mm 6 Qty
Bone plate 2.7mm 7 Qty
Bone plate 3.5/4.0mm 10 Qty
Cortical screws self-tapping 1.5mm 40 Qty
Cortical screws self-tapping 2.0mm 50 Qty
Cortical screws self-tapping 2.7mm 50 Qty
Cortical screws self-tapping 3.5mm 48 Qty
Cortical screws self-tapping 4.0mm-cancellous 32 Qty

Only for Veterinary use in USA.



COMPLETE ORTHOPEDIC SYSTEM FOR SMALL AUTOCLAVES



This is a complete kit of 1.5/2.0/2.7/3.5/4.0mm orthopedic instruments, bone screws, and bone plates for small autoclaves. Everything is constructed of high quality stainless steel. The cases are made of strong aluminium. All screws are self-tapping with a hex head.

1 Aluminum Case (Approx. Dimensions: 11 3/4"L x 7"W x 4"H)

1 Screw Case 1.5/2.0/2.7/3.5/4.0mm (Approx. Dimensions: 9 1/2"L x 4 1/2"W x 2 5/8"H)

1.1 Drill Bit (2)

1.5 x 80mm Drill Bit (2)

2.0 x 80mm Drill Bit (2)

2.5 x 115mm Drill Bit (2)

3.2 x 150mm Drill Bit (2)

1.1mm Neutral & Loaded Drill Guide

1.5mm Neutral & Loaded Drill Guide

2.0mm Neutral & Loaded Drill Guide

2.5mm Neutral & Loaded Drill Guide

3.2mm Neutral & Load Drill Guide

1.1mm & 1.5mm Tap Sleeve

2.5mm & 3.5mm Tap Sleeve

3.2mm & 4.0mm Tap Sleeve

2.0mm & 2.7mm Tap Sleeve

HA 1.5mm Tap

HA 2.0mm Tap

HA 2.7mm Tap

HA 3.5mm Tap

HA 4.0mm Tap

Small Screw Forceps

Hex Head Screw Driver SW2.5mm

Hex head Screwdriver 1.5mm

Large 70mm Depth Gauge

Small 40mm Depth Gauge

Large 6.0mm Counter Sink

Small 4.5mm Counter Sink

Small Plate benders (pair)

Large Plate benders (pair)

Bone plate 1.5mm (5)

Bone plate 2.0mm (6)

Bone plate 2.7mm (7)

Bone plate 3.5/4.0mm (10)

Cortical screws self-tapping hex head 1.5mm (40)

Cortical screws self-tapping hex head 2.0mm (50)

Cortical screws self-tapping hex head 2.7mm (50)

Cortical screws self-tapping hex head 3.5mm (48)

Cancellous screws self-tapping hex head 4.0mm (27)

For Veterinary Use Only in the USA.



COMPRESSION BONE PLATING KIT 3.5/4.0MM




Fantastic quality and very well packed! Excellent quality orthopedic implants/instruments is made of quality stainless steel. All the screws are self tapping-hex head. This kit can be only for veterinary use in the USA.

This kit can be adjusted according to customer's requirements. Instruments and implants included in this Compression Bone Plating Kit are listed below.

- Drill bit 3.5 2 Qty
- Drill bit 2.5 2 Qty
- HB4.0 Tap Sleeve & 3.5 Drill Guide
- HA3.5 Tap
- HB4.0 Tap
- Depth Gauge
- Hex Head Screw Driver SW2.5
- Screw Holding Sleeve
- Screw Holding Forceps
- Neutral & Loaded Drill Guide /2.5 HA3.5
- HA3.5 Tap Sleeve & 2.5 Drill Guide
- Plate Benders
- Periosteal Elevator (Round Head)
- Periosteal Elevator (Flat Head)
- Pry Bar
- Reduction Forceps with point
- Reduction Forceps
- Self-centering forceps 2 Qty
- Broken Screw Exteriorizer
- Countersink Drill
- Case (with screw rack)
- Compression Bone plates 3.5 mm 5 Qty
- Cortical screws 3.5 mm 32 Qty
- Cancellous screws 4mm 18 Qty

Only for Veterinary use in USA.

CUSTOM MADE ORTHOPEDIC KIT



Custom Designed Veterinary Orthopedic System

We can provide a custom made system according to your needs. We offer all kinds of orthopedic instruments, implants and instrument boxes that will fit different kinds of autoclaves. Implants are made of stainless steel or titanium. Bellow are your choices.

Drill Bits, all Jacob Chuck

Size 1.0mm, 1.3mm, 1.5mm, 1.7mm, 2.0mm, 2.5mm, 2.7mm, 3.2mm, 3.5mm, 4.0mm

Bone Screws, all hex head, cortical, self tapping

1.5mm (6-20mm)

2.0mm (8-30mm)

2.7mm (10-30mm)

3.5mm (10-50mm)

Bone Screws, cancellous, fully threaded

4.0mm (12-50mm)

6.5mm (10-30mm)

Screw Drivers, all sizes

quick coupling

fiber handle

Bone Plates

1.5mm (3-12 hole length)

2.0mm (3-12 hole length)

2.7mm (3-12 hole length)

3.5mm (3-12 hole length)

Bone Plates

4.5mm (5-20 hole length)

6.5mm (5-20 hole length)

Neutral and load drill guides, all size

Five different sizes of instrument cases, according you need



INSTRUMENT KIT 4.5/6.5MM



Instruments only

- Neutral & Loaded drill guide 3.2
- Depth gauge 30
- Countersink Drill 8
- Broken screw exteriorizer
- Plate bending forceps 2x1
- Periosteal elevator Flat 16
- Periosteal elevator Round 8
- Screwdriver (hex) + driver nut sw 3.5
- Hohmann retractor 2x1
- HA 4.5 Tap
- HA 6.5 Tap
- Skid exteriorizer sw 3.5
- Reduction forceps
- Curved reduction forceps
- Self - centering forceps (two)
- Drill bit - 3.2x150 (two)
- Drill bit - 4.5x150 (two)
- C-Shaped Guide 3x1
- Drill Bushing 3.2/4.5
- Drill bushing 4.5/6.5
- Screw holding forceps
- Case with 4.5/6.5 screw rack

Veterinary Use only in USA.

INSTRUMENT SYSTEM-KIT 3.5MM/4.0MM



Instrument System Kit Orthopedic

The metal case is aluminum. It is very good quality and very strong. It comes with a cover. The quality is excellent and the manufacturer has CE and other certificates. All screws are self tapping, hex head with standard size of the head. Packed in nice aluminum. The dimension is 15.5" x 9.5" x 4" (L x W x H) and the case has 3 levels for storing the instruments. Below are list of instruments included in the Instrument System-Kit.

Drill bit 3.2 = 2

Drill bit 2.5 = 2

Neutral Loaded Drill Guide 2.5 HA3.5 = 1

HA3.5 Tap Sleeve 2.5 Drill Guide = 1

HB4.0 Tap Sleeve 3.0 Drill Guide = 1

HA3.5 Tap = 1

HB4.0 Tap = 1

Depth Gauge = 1

Hex Head Screw Driver SW2.5 = 1

Screw Holding Sleeve = 1

Plate Benders = 1

Periosteal Elevator (Round Head) = 1

Periosteal Elevator (Flat Head) = 1

Pry Bar = 1

Reduction Forceps with point = 1

Reduction Forceps = 1

Self-centering forceps = 2


Countersink Drill = 1

Case (screw rack) 3.5/4.0mm = 1

For veterinary USA only in USA



LARGE FRAGMENTS KIT 4.5/6.5MM



Veterinary Orthopedic 4.5 - 6.5mm System.

This excellent quality system offers plates, screws and instruments. All screws are self tapping, hex head with standard size of the head. Packed in nice aluminum. This kit can be only for veterinary use in the USA.

Below are list of (plates, screws, instruments) included in the Veterinary Orthopedic 4.5 - 6.5mm System.

Neutral & Loaded drill guide 3.2

Depth gauge 30

Countersink Drill 8

Broken screw exteriorizer

Plate bending forceps 2x1

Periosteal elevator Flat 16

Periosteal elevator Round 8

Screwdriver (hex) + driver nut sw 3.5

Hohmann retractor 2x1

HA 4.5 Tap

HA 6.5 Tap

Skid exteriorizer sw 3.5

Reduction forceps

Curved reduction forceps

Self - centering forceps (two)

Drill bit - 3.2x150 (two)

Drill bit - 4.5x150 (two)

C-Shaped Guider 3x1

Drill Bushing 3.2/4.5

Drill bushing 4.5/6.5

Screw holding forceps

Case with 4.5/6.5 screw rack

Plates - different sizes for 4.5 and 6.5 screws = 10 (5-14 holes 92-254mm)

Bone screws HA 4.5 (cortical) and 6.5(cancellous) = 90



MEDIUM KIT 2.7/3.5/4.0MM VETERINARY



High quality orthopedic implants/instruments. The retail price of this pack is more than 2000 USD. You can get it here for a good price and the quality is excellent. All the instruments and implants are made from the same material used for human medicine. This set has all the instruments as others do, plus this pack has 10 more instruments.

Excellent quality. The material is all stainless steel. The metal case is aluminum and it comes with a very strong cover. Dimensions is 15.5" x 9.5" x 4" (L x W x H). Also includes a screw rack. As you can see from the picture this is packed with instruments and screws. The 3.5 mm plates are very strong and individually packed. Optional-Titanium implants.

If you need to know more about the sizes of the screws and plates, please contact us. We are sure that you will like it. If you are not satisfied, you can return it and we will be glad to refund you.

System Includes: Drill Bits

(2) 2.0mm

(2) 2.5mm

(2) 3.2mm

Neutral & Load Drill Guides

(1) 2.0mm

(1) 2.5mm

(1) 3.2mm

Combined Drill Guides & Tap Sleeves

(1) 2.0mm/2.7mm

(1) 2.5mm/3.5mm

(1) 3.2mm/4.0mm

Bone Taps

(1) 2.7mm

(1) 3.5mm

(1) 4.0mm

Other Instruments

(1) Depth Gauge

(1) 2.5mm Hex Head Screwdriver with Screw Holding Sleeve

(1) Small Screw Holding Forceps

(1) Pair of Large Plate Benders

(1) Pair of Small Plate Benders

(1) Periosteal Elevator (Round Head)

(1) Periosteal Elevator (Flat Head)

(1) Broken Screw Extractor

(1) Reduction Forceps

(1) Reduction Forceps with Point

(2) Self-Centering Forceps

(1) Small Countersink

(1) Large Countersink

Bone Plates

(4) 2.7mm

(9) 3.5/4.0mm

Bone Screws

(30) 2.7mm Cortical Self-Tapping Hex Head


(60) 3.5mm Cortical Self-Tapping Hex Head

(18) 4.0mm Cancellous Self-Tapping Hex Head

For veterinary use only in USA.



MICRO KIT 1.5/2.0/2.7MM



Micro Kit 1.5/2.0/2.7mm Orthopedic kit


Customers can request different sizes and numbers of plates and screws if needed. Excellent in quality. It is made of all stainless steel material. All the instruments and implants are made from the same material used for human medicine.

Instruments and implants included in this complete kit are listed below.

- Bone drill bit 1.1 x 80mm = 5 pcs
- Bone drill bit 1.5 x 80mm = 5 pcs
- Bone drill bit 2.0 x 80mm = 5 pcs
- Tap 1.5mm = 1 pcs
- Tap 2.0mm = 1 pcs
- Tap 2.7mm = 1 pcs
- Screwdriver 1.5/2.0mm = 1 pcs
- Screw holding sleeve 1.5mm = 1 pcs
- Screwdriver 2.7mm = 1 pcs
- Screw holding sleeve 2.7mm = 1 pcs
- Guider 1.1/1.5mm = 1 pcs
- Guider 1.5/2.0mm = 1 pcs
- Guider 2.0/2.7mm = 1 pcs
- Counter drill 4.0mm = 1 pcs
- Straight quick-coupling handle 4.5mm = 1 pcs
- T-shape quick-coupling handle 4.5mm = 1 pcs
- Screwdriver with handle 2.6mm = 1 pcs
- Depth detector 40mm = 1 pcs
- Hook = 2 pcs
- De-coupling = 1 pcs
- Periosteal dissector 3mm = 1 pcs
- Periosteal dissector 5mm = 1 pcs
- Screw holder = 1 pcs
- Plate forcep = 1 pcs
- Plate bender 2*1 = 1 set
- Elevator 3*1 = 1 pcs
- Bent reduction forcep = 1 pcs
- Sharp reduction forcep = 1 pcs
- Plate nipper = 1 pcs
- Multi-functional forcep = 1 pcs
- Bone plate box 1.5mm = 1 set
- Bone plate box 2.0mm = 1 set
- Bone plate box 2.7mm = 1 set
- Case with lid 15" x 9" x 3"
- Bone screws 1.5mm self tapping = 40 screws (6-30mm)
- Bone screws 2.0mm self tapping = 55 screws (6-30mm)
- Bone screws 2.7mm self tapping = 60 screws (8-30mm)
- Bone plate 1.5mm = 7
- Bone plate 2.0mm = 10
- Bone plate 2.7mm = 10

USA can be for Veterinary use only.


MICRO-KIT ONLY INSTRUMENTS



Micro-kit only instruments offers superior in quality, style and design. These instruments are easy to use. Great for your orthopedic surgical instrument needs. Below are the instruments.

Bone drill bit 1.1 x 80mm = 5 pcs
 Bone drill bit 1.5 x 80mm = 5 pcs
 Bone drill bit 2.0 x 80mm = 5 pcs
 Tap 1.5mm = 1 pcs
 Tap 2.0mm = 1 pcs
 Tap 2.7mm = 1 pcs
 Screwdriver 1.5/2.0mm = 1 pcs
 Screw holding sleeve 1.5mm = 1 pcs
 Screwdriver 2.7mm = 1 pcs
 Screw holding sleeve 2.7mm = 1 pcs
 Guider 1.1/1.5mm = 1 pcs
 Guider 1.5/2.0mm = 1 pcs
 Guider 2.0/2.7mm = 1 pcs
 Counter drill 4.0mm = 1 pcs
 Straight quick-coupling handle 4.5mm = 1 pcs
 T-shape quick-coupling handle 4.5mm = 1 pcs
 Screwdriver with handle 2.6mm = 1 pcs
 Depth detector 40mm = 1 pcs
 Hook = 2 pcs
 De-coupling = 1 pcs
 Periosteal dissector 3mm = 1 pcs
 Periosteal dissector 5mm = 1 pcs
 Screw holder = 1 pcs
 Plate forcep = 1 pcs
 Plate bender 2*1 = 1 set
 Elevator 3*1 = 1 pcs
 Bent reduction forcep = 1 pcs
 Sharp reduction forcep = 1 pcs
 Plate nipper = 1 pcs
 Multi-functional forcep = 1 pcs
 Bone plate box 1.5mm = 1 set
 Bone plate box 2.0mm = 1 set
 Bone plate box 2.7mm = 1 set
 Case with lid 15" x 9" x 3"
 In USA can be for Veterinary use only.

MINI KIT 1.5/2.0/2.7MM




This Orthopedic System 1.5/2.0/2.7mm includes quality instruments and implants. We provide different sizes and numbers of plates and screws with the request of the customer's needs.

Comes with plates, screws and instruments. All screws are self-tapping. Instruments and implants included in this Mini Kit 1.5/2.0/2.7mm kit are listed below.

Drill bits 1.3/1.7/2.0 = 6
 Guider HA 2.0 = 1
 Tap HA 2.0 = 1
 Countersink drill HA 2.0 = 1
 Periosteal elevator 3 mm = 1
 Screw holding forceps = 1
 Tip pull hook = 1
 Periosteal elevator 5mm = 1
 Skin retractor = 1
 Hofmann retractor 3x1 = 1
 Depth Gauge = 1
 Plate Bending wrench = 1
 Bone Bending forceps pair = 1
 Plate holder = 1
 Reduction forceps = 1
 Curred reduction forceps = 1
 Bone Rogers = 1
 Wire cutter = 1
 Screw driver = 1
 Metal case with cover = 1
 Bone screws 2.0 mm = 50
 Bone screws 2.7 mm = 50
 Bone plate 2mm = 7
 Bone plate 2.7mm = 7
 Bone plate 1.5mm = 7
 Bone screws 1.5mm = 40

Only for Veterinary use in USA.

STARTER KIT 2.7/3.5MM



Veterinary orthopedic Starter Kit 2.7/3.5mm

Excellent quality at an affordable price!

Orthopedic system for animal surgeries. This set has all the instruments for a starter kit. The material is all stainless steel. It is made from the same material used for human medicine.

Instruments and implants included in this Starter Kit 2.7/3.5mm are listed below.



Drill bit 2.0 = 2
 Drill bit 2.5 = 2
 Drill bit 2.7 = 2
 Neutral and load Drill guide 2.0mm = 1
 Neutral-Loaded Drill Guide HA 2.5mm = 1
 HA3.2 Tap Sleeve 2.5 Drill Guide = 1
 HB4.0 Tap Sleeve 3.5 Drill Guide = 1
 HA3.5 Tap = 1
 HB4.0 Tap = 1
 Depth Gauge = 1
 Hex Head Screw Driver SW2.5mm = 1
 Screw Holding Sleeve = 1
 Plate Benders pair-modern design = 1
 Periosteal Elevator(Round Head) = 1
 Periosteal Elevator (Flat Head) = 1
 Reduction Forceps = 1
 Self-centering forceps = 1
 Countersink Drill = 1
 Case with screw rack - modern
 Bone plate 2.7mm = 5
 Bone plate 3.5mm = 7
 Cortical screws 2.7mm = 30
 Cortical screws 3.5 mm = 64
 (Screws self-tapping, hex head)

For veterinary USA only in USA