

Otyker orthopaedics

Otyker Orthopaedics Technology Co. Ltd.

Phone: +852 8191 3424

Email: contact@otyker.co

Skype: otykerco

Address: 212Rm, No.592 Tong He Road,
Qingdao 266035 China



BONE PLATE 1.5MM T SHAPE

T Shape bone plate 1.5mm. Made of stainless steel. Light in weight, but very strong. Sizes available are listed below. 3 holes 4 holes 5 holes 6 holes 7 holes 8 holes 9 holes



BONE PLATE 2.0MM T SHAPE

T Shape bone plate for 2.0mm. Very strong. Made of stainless steel. Sizes available are listed below. 3 holes 4 holes 5 holes 6 holes 7 holes 8 holes 9 holes



BONE PLATE 2.7MM BY GPC

Excellent quality, very strong. Made of stainless steel material. It is FDA approved. Sizes offers are list below. 3 holes 4 holes 5 holes 6 holes 7 holes 8 holes 9 holes 10 holes 11 holes 12 holes 13 holes



BONE PLATE 2.7MM T SHAPE

Very Strong - Bone plate for 2.7mm T shape. Made from stainless steel. Sizes available are listed below. 3 holes 4 holes 5 holes 6 holes 7 holes 8 holes 9 holes 10 holes 11 holes 12 holes



BONE PLATE 3.5MM BY GPC

Excellent quality, very strong. Made of stainless steel material. It is FDA approved. Sizes offers are list below. 3 holes 4 holes 5 holes 6 holes 7 holes 8 holes 9 holes 10 holes 11 holes 12 holes




BONE PLATE 3.5MM T SHAPE

Bone plate for 3.5mm T shape. Made of stainless steel, it is very strong. Contact us for size.




BONE PLATE 4.5/6.5MM

Straight plates - Bone plate 4.5/6.5mm. Made of stainless steel, it is very strong. Will fit to 4.5mm and 6.5mm screws. Sizes offers are list below. 4 holes 5 holes 6 holes 7 holes 8 holes 9 holes 10 holes 11 holes 12 holes 13 holes 14 holes 15 holes



BONE PLATE 6.5 MM PREMIUM QUALITY

Veterinary orthopaedic bone plate 6.5 mm. Excellent quality and very strong. It is FDA and CE approved. Made from stainless steel. Has great customer reviews! (Note: Picture is of size 3.5mm).




BONE SCREW 4.5MM

Cancellous fully threaded. Made of stainless steel, it is very strong. Has great customer reviews.




BONE SCREW 6.5MM

Cancellous fully threaded. Hex head. Made of stainless steel, very strong. Has excellent customer reviews. Sizes available are listed below. 30mm 35mm 40mm 45mm 50mm 55mm 60mm 65mm 70mm 75mm 80mm 85mm 90mm




BONE SCREW CANCELLOUS 3.5MM

Hex head cancellous - 2.5mm head socket. Excellent quality. Made of stainless steel.




BONE SCREW CANCELLOUS 4.0MM

Bone screw cancellous 4.0mm. Excellent quality, very strong. Made of stainless steel. Sizes available are listed below. 12mm 14mm 16mm 18mm 20mm 22mm 24mm 26mm 28mm 30mm 32mm 34mm 36mm 38mm 40mm 42mm 44mm 46mm 48mm 50mm



BONE SCREW CANCELLOUS 4.5MM

Bone screw cancellous 4.5mm. Excellent quality, very strong. Made of stainless steel.



BONE SCREW CANCELLOUS 6.0 MM

Bone screw cancellous 6.0mm. Excellent quality, very strong. Made of stainless steel.



ORTHOPAEDIC BONE PLATES 1.5MM PREMIUM QUALITY

Excellent quality - veterinary orthopaedic implants. Animal bone plate 1.5mm. Very strong. Made of stainless steel. Has excellent customer reviews! Sizes available are listed below. 4 holes 5 holes 6 holes 7 holes 8 holes 9 holes 10 holes 11 holes 12 holes 13 holes



ORTHOPAEDIC BONE PLATES 1.5MM SMALL

Veterinary orthopaedic implants - Animal bone plate for 1.5mm. Excellent quality and very strong. Made of stainless steel. Very thin(1mm) plates and very good for small animals and birds. Easy to band even with fingers. Sizes available are listed below. 3 holes 4 holes 5 holes 6 holes 7 holes 8 holes 9 holes 10 holes 11 holes 12 holes

ORTHOPAEDIC BONE PLATES 2.0MM PREMIUM QUALITY

Special veterinary plates designed for small animals. It is for 2.0mm screws. Very strong and made from high quality medical grade stainless steel (1.5mm thick). Not sterile. Sizes offers are list below. 4 holes 5 holes 6 holes 7 holes 8 holes 9 holes 10 holes 11 holes 12 holes 13 holes

ORTHOPAEDIC BONE PLATES 2.0MM SMALL

Animal bone plate for 2.0mm. Excellent veterinary orthopaedic implants-very strong. Made of stainless steel. These plates are thin (1mm) and extremely good for small animals and birds. Very easy to band. Sizes available are listed below. 3 holes 4 holes 5 holes 6 holes 7 holes 8 holes 9 holes 10 holes 11 holes 12 holes 13 holes

ORTHOPAEDIC BONE PLATES 2.7MM PREMIUM QUALITY

Veterinary orthopaedic animal implants for 2.7mm. Excellent quality, very strong. Made of stainless steel. Thickness is 2.5mm. Holes is oval. Has excellent customer reviews! Sizes offers are list below. 4 holes 5 holes 6 holes 7 holes 8 holes 9 holes 10 holes 11 holes 12 holes 13 holes

ORTHOPAEDIC BONE PLATES 3.5MM PREMIUM QUALITY

Veterinary orthopaedic bone plate 3.5 mm. Excellent quality and very strong. It is FDA and CE approved. Made from stainless steel. Screw size: 3.5mm cortical, 4.0mm cancellous. Sizes available are listed below. 3 holes 4 holes 5 holes 6 holes 7 holes 8 holes 9 holes 10 holes 11 holes 12 holes 13 holes 14 holes 15 holes 16 holes

ORTHOPAEDIC BONE PLATES 4.5/6.0MM

Veterinary Orthopaedic LD-DCP bone plate for 4.5/6.0 mm screws. Made from stainless steel. Excellent quality and very strong. Has great customer reviews! Sizes available are listed below. 3 holes 4 holes 5 holes 6 holes 7 holes 8 holes 9 holes 10 holes 11 holes 12 holes 13 holes

ORTHOPAEDIC BONE SCREWS 1.5MM SELF TAPPING HEX HEAD

Premium quality 1.5mm orthopaedic bone screws. Hex head fully threaded. Made of stainless steel. Self Tapping. Has excellent customer review! Sizes available are listed below. 6mm 8mm 10mm 12mm 14mm 16mm 18mm 20mm

ORTHOPAEDIC BONE SCREWS 2.0MM SELF TAPPING HEX HEAD

2.0mm orthopaedic bone screws. Very good quality. Self tapping. Cortical hex head fully threaded. Made of stainless steel. Sizes available are listed below. 6mm 8mm 10mm 12mm 14mm 16mm 18mm 20mm 22mm 24mm 26mm 28mm 30mm

ORTHOPAEDIC BONE SCREWS 2.7MM SELF TAPPING HEX HEAD

Orthopaedic Bone screws 2.7mm- Hex head fully threaded. Made of stainless steel. Excellent quality. Self Tapping. Sizes available are listed below. 10mm 12mm 14mm 16mm 18mm 20mm 22mm 24mm 26mm 28mm 30mm



from
\$5

ORTHOPAEDIC BONE SCREWS 3.5MM SELF TAPPING HEX HEAD

Quality veterinary Orthopaedic screws 3.5mm is made of stainless steel. Excellent quality. Hex head fully threaded. Self Tapping. Sizes available are listed below. 10mm 12mm 14mm 16mm 18mm 20mm 22mm 24mm 26mm 28mm 30mm 32mm 34mm 36mm 38mm 42mm 44mm 46mm 48mm 50mm

from
\$99

RECONSTRUCTION BONE PLATE

Reconstruction Bone Plate 15hl - Made of stainless steel. Very strong. Comes in sizes of 2.7mm or 3.5mm. 2.7mm is about 5 inches long 3.5mm is 7 inches long

from
\$99

TPLO BONE PLATE 2.7MM

Premium quality - TPLO Bone plate for 2.7 screws. Made from stainless steel. Very strong. Has great customer reviews.

from
\$149

VETERINARY ORTHOPAEDIC TPLO PLATE 3.5MM

Quality TPLO plates for animal surgeries. For screws 3.5mm. Made of stainless steel material. Comes with 3 holes on the body or with 4 holes.